



东莞中力仪器科技有限公司

ZhongLi Instrument Technology Co., Ltd.

Introduction

Name: ZL-1026 High Temperature Muffle Lab Furnace Manufacturer



Tel: [+86-769-33657986](tel:+86-769-33657986)

Fax: [+86-769-33657987](tel:+86-769-33657987)

Website: www.zltester.com

E-mail: Zhongli@zltester.com

Add: No.62, Xinhe Management of Area, Wanjiang District, Dongguan, Guangdong, China



Specifications

- 1.High Temperature Muffle Lab Furnace
- 2.Power rating:5Kw
- 3.temperature:1200C (can be customized)
- 4.oven size:200*300*120mm (can be customized)
- 5.3 years warranty

Product Description:

High temperature muffle lab furnace used on analysis element and measuring the function of heating and heating treatment of General small steel hardening, annealing, tempering. High - temperature furnace can also be used for metal, ceramic sintering, dissolving, analysis of high temperature heating

Main Feature:

The lab muffle furnace with modern high - tech technology, temperature error ± 2 °C, precision temperature control ± 2 °C. Main parts adopt imported material, shell is advanced CNC machine tooling, imported luxury, beautiful two - color epoxy powder electrostatic painting processes, furnaces using a special light material and it. The indicators have reached the international advanced level.

Main Usage:

Primarily for ceramics, glass, enamel, and inorganic pigments, industrial and mining enterprises and scientific research on abrasive fritz glaze, glass, low - temperature flux, used in enamel and bond test, also for alloys, composite materials, melting, physical and chemical experiments.

Apply Industry:

High temperature muffle lab furnace can be use on determination, metallographic laboratory for analysis, research, or industrial and mining enterprises for high - temperature sintering and metal work in the air for normalizing, annealing, hardening, for the Laboratory for universities, scientific research and environmental protection unit.

Technical parameters:

Model	ZL-1026 High temperature muffle lab furnace
Electric power	220V

Tel: [+86-769-33657986](tel:+86-769-33657986)

Fax: [+86-769-33657987](tel:+86-769-33657987)

Website: www.zltester.com

E-mail: Zhongli@zltester.com

Add: No.62, Xinhe Management of Area, Wanjiang District, Dongguan, Guangdong, China



东莞中力仪器科技有限公司

ZhongLi Instrument Technology Co., Ltd.

Power rating	5Kw
Power temperature	1200°C
Size of working areas	200*300*120mm
Around outer size	480*460*630mm
Heating elements	silicon carbon stick
Temperature control	thermonuclear(WREU)
Support table	Theorist surge console

Tel: [+86-769-33657986](tel:+86-769-33657986)

Fax: [+86-769-33657987](tel:+86-769-33657987)

Website: www.zltester.com

E-mail: Zhongli@zltester.com

Add: No.62, Xinhe Management of Area, Wanjiang District, Dongguan, Guangdong, China
